## anti-human CD62L PE-conjugated

PE - conjugated monoclonal antibody LT-TD180 to human CD62L

Cat-No: **21279624** 500 μl

Clone: LT-TD180

**Specificity:** The antibody LT-TD180 reacts with CD62L (L-selectin), a 74-95 kDa single chain type I glycoprotein expressed on most peripheral blood B lymphocytes, T lymphocytes, monocytes and granulocytes; it is also present on a subset of NK cells and certain hematopoietic malignant cells.

Isotype subclass: Mouse IgG1

**Form:** The purified antibody is conjugated with R-Phycoerythrin (PE) under optimum conditions. The reagent is adjusted for direct use. No reconstitution is necessary.

**Expiration date:** The reagent is stable until the expiry date stated on the vial label

Physical state: Liquid

Buffer/Additives/Preservative: PBS containing 1 % BSA and 0.09 % sodium azide (pH 7.2)

Storage conditions: Store at 4 °C. Do not freeze. Avoid prolonged exposure to light.

**Application:** Flow Cytometry

## References:

\*Von Andrian UH and others: 1992 Oct;263(4 Pt 2):H1034-44.

- \*von Andrian UH and others: 1993 Jul 1;82(1):182-91.
- · \*Simon SI and other: 1995 Aug 1;155(3):1502-14.
- \* \*Ding Z and other: 2003 Jun 1;101(11):4245-52.
- \*Ramachandran V and other: 2004 Sep 14;101(37):13519-24.

**Background:** CD62L (L-selectin) is an adhesion glycoprotein that is constitutively expressed on the cell surface of leukocytes and mediates their homing to inflammatory sites and peripheral lymph nodes by enabling rolling along the venular wall. CD62L is also involved in activation-induced neutrophil aggregation. Activation-dependent CD62L shedding, however, counteracts neutrophil rolling. CD62L has also signaling roles including enhance of chemokine receptor expression. Similarly to CD62P, the major ligand of CD62L is PSGL-1 (P-selectin glycoprotein ligand-1).

**Warning:** Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32). Azide compounds should be flushed with large volumes of water during disposal to avoid deposits in lead or copper plumbing where explosive conditions can develop.

This material is offered for <u>research use only</u>. Not for use in human. For in vitro use only. ImmunoTools will not be held responsible for patent infringement or other violations that may occur with the use of our products.

## ImmunoTools Excellent Quality - Advantageously priced

Gladiolenweg 2; 26169 Friesoythe; Germany

phone:+49-(0)4491-400997, fax:+49-(0)4491-400998, info@immunotools.com

www.immunotools.com