

Mouse IgG1 Isotype control PerCP-conjugated

PerCP-conjugated monoclonal antibody 1F8 to IgG1

Cat-No: **21815015**

500 µl

Clone: 1F8

Specificity: This reagent can be used as a mouse IgG1 isotype control. It is obtained from a fused hybridoma by spleen cell from Balb/C mouse immunized by tylosin and KLH, and myeloma SP2/0 from Balb/C mouse. The antibody is specific for tylosin, may bind non-specifically to KLH. This hapten is not expressed on human cells or human cell-lines.

Isotype subclass: Mouse IgG1

Form: The purified antibody is conjugated with Peridinin Chlorophyll (PerCP) under optimum conditions. The conjugate is purified by size-exclusion chromatography and adjusted for direct use. No reconstitution is necessary.

Physical state: Liquid

Buffer/Additives/Preservative: PBS containing 1% BSA and 0.09% sodium azide (pH 7.2).

Expiration date: The reagent is stable until the expiry date stated on the vial label.

Storage conditions: Store at 4 °C. Avoid prolonged exposure to light.

Application: Flow Cytometry

Background: The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.

Warning: Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32). Azide compounds should be flushed with large volumes of water during disposal to avoid deposits in lead or copper plumbing where explosive conditions can develop.

This material is offered for **research use only**. Not for use in human. For in vitro use only. ImmunoTools will not be held responsible for patent infringement or other violations that may occur with the use of our products.

ImmunoTools Excellent Quality - Advantageously priced

Gladiolenweg 2; 26169 Friesoythe; Germany
phone: +49-(0)4491-400997, fax: +49-(0)4491-400998, info@immunotools.com
www.immunotools.com