

## Mouse IgG1 control APC-conjugated

**APC-** conjugated monoclonal antibody PPV-06 for negative control

Cat-No: **21275516**

500 µl

**Clone:** PPV-06

**Specificity:** This mouse IgG1 monoclonal antibody (clone PPV-06) reacts with undefined epitope on a plant pathogen.

**Negativ Species:** Human, Porcine, Mouse and Rat

**Isotype subclass:** Mouse IgG1

**Form:** The purified antibody is conjugated with cross-linked Allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

**Physical state:** Liquid

**Buffer/Additives/Preservative:** PBS containing BSA and 15mM sodium azide (pH 7.4)

**Expiration date:** The reagent is stable until the expiry date stated on the vial label.

**Storage conditions:** Store at 4 °C. Avoid prolonged exposure to light.

**Application:**

Flow Cytometry  
Control Experiments

**Background:** The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.

**Warning:** Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32). Azide compounds should be flushed with large volumes of water during disposal to avoid deposits in lead or copper plumbing where explosive conditions can develop.

This material is offered for **research only**. Not for use in human. For in vitro use only. ImmunoTools will not be held responsible for patent infringement or other violations that may occur with the use of our products.

**ImmunoTools** Excellent Quality - Advantageously priced

Altenoyther Str. 10; 26169 Friesoythe; Germany  
phone:+49-(0)4491-400997, fax:+49-(0)4491-400998, info@immunotools.com  
www.immunotools.com