

anti-human Anti-Fitc

Monoclonal Antibody F4/1 to Anti-Fitc (Human)

Cat-No: 21153331

100 µg in 100 µl

Clone: F4/1

Specificity: The Anti-FITC (fluorescein isothiocyanate) monoclonal antibody specifically detects both the free and protein conjugated (either soluble or cell-bound) form of fluorescein or carboxyfluorescein.

Isotype subclass: Mouse IgG1_κ

Form: Purified

Physical state: Liquid

Buffer/Additives/Preservative: PBS containing 0.09 % sodium azide (pH 7.4).

Expiration date: The reagent is stable until the expiry date stated on the vial label.

Storage conditions: Store at 4°C. For long term storage, aliquot and freeze unused portions at -20°C. Avoid repeated freeze/thaw cycles.

Application: This antibody is suitable for use in immunohistochemistry (acetone-fixed frozen section), flow cytometry and immunoprecipitation of cells or molecules labeled with fluorescein derivatives (FITC or CFSE). Also, this mAb is able to quench FITC fluorescence. There is no endogenous tissue staining with this antibody.

References: *Balogh P, Szekeres G, Nemeth P. 1994 Hapten-mediated identification of cell membrane antigens using an anti-FITC monoclonal antibody. J Immunol Methods. 1994 169:35-40

Warning: Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32). Azide compounds should be flushed with large volumes of water during disposal to avoid deposits in lead or copper plumbing where explosive conditions can develop.

This material is offered for research only. Not for use in human. For in vitro use only. ImmunoTools will not be held responsible for patent infringement or other violations that may occur with the use of our products.

ImmunoTools Excellent Quality - Advantageously priced

Altenoyther Str. 10; 26169 Friesoythe; Germany
phone:+49-(0)4491-400997, fax:+49-(0)4491-400998, info@immunotools.com
www.immunotools.com